

ENDOCUFF VISION™ Compatibility Chart*

FUJI

Scope Model	ENDOCUFF VISION™ Model
EC-760R-V/L Adult	ARV110 (Blue)
EC-760ZP-VM/VL Adult	ARV110 (Blue)
EC-600WM/WI/WL New Series	ARV110 (Blue)
EC-590WM4/WI4/WL4/WI/WLY Adult	ARV110 (Blue)
EC-590ZW3 M/L Adult Zoom	ARV120 (Green)
EC-530FI/FL/FM Adult	ARV120 (Green)
EC-530WM3/WI3/WL3/WI	ARV120 (Green)
EC-450 Adult Series	ARV120 (Green)
EC-200MR Adult	ARV120 (Green)
EC-250 Adult Series	ARV120 (Green)
EC-201WI Adult	ARV120 (Green)
EC-530 MT/IT/LT Adult	ARV120 (Green)
EC-530 DM/DL Adult	ARV120 (Green)
EC-600 ZW/M/L Adult	ARV120 (Green)
EC-530MP/LP Pediatric	ARV130 (Purple)
EC-530LS Intermediate	ARV130 (Purple)
EC-450 Pediatric Series	ARV130 (Purple)
EC-450LP/LS5 Intermediate	ARV130 (Purple)
EC-450 Intermediate Series	ARV130 (Purple)

PENTAX

EC38-i10M/M2/F/F2/LNF Adult	ARV110 (Blue)
EC-3890TMK/MK2/FK/FK2/LK Adult	ARV110 (Blue)
EC-3890MZI/FZI/LZI/SI2/TFI Adult	ARV110 (Blue)
EC-3890FI/LI/FI2(HD)/Mi2 Adult	ARV110 (Blue)
EC-3890TLK Dual Channel Adult	ARV110 (Blue)
EC-3970LK Adult	ARV120 (Green)
EC-3880FK2 Adult	ARV120 (Green)
EC-3870LK/LZK/TLK/FK2/MK2 Adult	ARV120 (Green)
EC-3840FK2 Adult	ARV120 (Green)
EC-3830LK/P/FK2 Adult	ARV120 (Green)
EC-380LKP/FKP/FK2P/MK2PAdult	ARV120 (Green)
EC-38LV Adult	ARV120 (Green)
EC-3490FK Pediatric HD	ARV130 (Purple)
EC-3490 FI Pediatric HD	ARV130 (Purple)
EC-3490 MK/FK/LK/LI/TLI/TFI Pediatric	ARV130 (Purple)
EC-3470LK Pediatric	ARV130 (Purple)
EC-34-i10M/F/L Pediatric HD	ARV130 (Purple)
EC-3490TLF/MI	Not compatible
EC-2990LI	Not compatible

STORTZ

13924-PKS Adult	ARV120 (Green)
13925-PKS Adult	ARV120 (Green)
13907-PKS Pediatric	ARV130 (Purple)
13906-PKS Pediatric	ARV130 (Purple)
13905-PKS Adult	ARV140 (Orange)
13904-PKS Adult	ARV140 (Orange)

* The table provides an indication of ENDOCUFF VISION™ compatibility with colonoscopes. For confirmation of specific models, contact an Olympus representative.

Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer. Olympus is a registered trademark of Olympus Corporation, Olympus America Inc., and/or their affiliates. ENDOCUFF VISION™ is a registered trademark of KeyMed (Medical & Industrial Equipment) Ltd.